

# Marathon-Ready™ cDNA Protocol-at-a-Glance (PT1156-2)

Marathon cDNA amplification is a fairly complex, multiday procedure. Please read the *User Manual* before using this abbreviated protocol, and refer to it often for interpretation of results during the course of your experiments. The Protocol-at-a-Glance is provided for your convenience, but is not intended for first-time users.

## VI. Rapid Amplification of cDNA Ends (RACE)

1. Prepare enough PCR master mix for all of the reactions plus one additional tube. The same master mix can be used for both 5'- and 3'-RACE reactions. For each 50- $\mu$ l reaction, mix the following reagents:

36 $\mu$ l	H <sub>2</sub> O
5 $\mu$ l	10X cDNA PCR Reaction Buffer
1 $\mu$ l	dNTP Mix (10 mM)
1 $\mu$ l	Advantage® 2 Polymerase Mix (50X)
43 $\mu$ l	Final volume

Mix well by vortexing (without introducing bubbles) and briefly spin the tube in a microcentrifuge.

2. **For 5'-RACE:** prepare PCR reactions as shown in Table II.  
**For 3'-RACE:** prepare PCR reactions as shown in Table III.  
Add the components in the order shown in 0.5-ml PCR tubes.

**TABLE II: SETTING UP 5'-RACE REACTIONS**

Component	Test Tube #: Description:	1 Experimental Sample	2 G3PDH Pos. Ctrl *	3 GSP 1 + 2 Pos. Ctrl	4 AP1 only Neg. Ctrl	5 GSP1 only Neg. Ctrl
Marathon-Ready cDNA		5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l
AP1 Primer (10 $\mu$ M)		1 $\mu$ l	1 $\mu$ l	---	1 $\mu$ l	---
GSP1 (antisense primer; 10 $\mu$ M)		1 $\mu$ l	---	1 $\mu$ l	---	1 $\mu$ l
GSP2 (sense primer; 10 $\mu$ M)		---	---	1 $\mu$ l	---	---
Control 5'-RACE G3PDH Primer (10 $\mu$ M)		---	1 $\mu$ l	---	---	---
H <sub>2</sub> O		---	---	---	1 $\mu$ l	1 $\mu$ l
Master Mix		43 $\mu$ l	43 $\mu$ l	43 $\mu$ l	43 $\mu$ l	43 $\mu$ l
Final volume		50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l

\* The G3PDH Positive Control should generate a 1.09-kb product.

**TABLE III: SETTING UP 3'-RACE REACTIONS**

Component	Test Tube #: Description:	1 Experimental Sample	2 G3PDH Pos. Ctrl †	3 GSP 1 + 2 Pos. Ctrl	4 AP1 only Neg. Ctrl	5 GSP2 only Neg. Ctrl
Marathon-Ready cDNA		5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l
AP1 Primer (10 $\mu$ M)		1 $\mu$ l	1 $\mu$ l	---	1 $\mu$ l	---
GSP2 (sense primer; 10 $\mu$ M)		1 $\mu$ l	---	1 $\mu$ l	---	1 $\mu$ l
GSP1 (antisense primer; 10 $\mu$ M)		---	---	1 $\mu$ l	---	---
Control 3'-RACE G3PDH Primer (10 $\mu$ M)		---	1 $\mu$ l	---	---	---
H <sub>2</sub> O		---	---	---	1 $\mu$ l	1 $\mu$ l
Master Mix		43 $\mu$ l	43 $\mu$ l	43 $\mu$ l	43 $\mu$ l	43 $\mu$ l
Final volume		50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l

† The G3PDH Positive Control should generate a 1.2-kb product.

3. Overlay the contents of each tube with 2 drops of mineral oil and place caps firmly on each tube.

**Note:** This is not necessary if you are using a hot-lid thermal cycler.

4. Commence thermal cycling using one of the following programs (programs 1 and 2 work with the Control G3PDH and AP1 Primers):

**Program 1** (preferred; use if GSP  $T_m > 70^\circ\text{C}$ ):

PE DNA Thermal Cycler 480:

- 94°C for 1 min
- 5 cycles:
  - 94°C 30 sec
  - 72°C 4 min<sup>‡</sup>
- 5 cycles:
  - 94°C 30 sec
  - 70°C 4 min<sup>‡</sup>
- 20–25 cycles:
  - 94°C 20 sec
  - 68°C 4 min<sup>‡</sup>

PE GeneAmp Systems 2400/9600  
(or hot-lid thermal cycler):

- 94°C for 30 sec
- 5 cycles:
  - 94°C 5 sec
  - 72°C 4 min<sup>‡</sup>
- 5 cycles:
  - 94°C 5 sec
  - 70°C 4 min<sup>‡</sup>
- 20–25 cycles:
  - 94°C 5 sec
  - 68°C 4 min<sup>‡</sup>

**Program 2** (if GSP  $T_m = 60\text{--}65^\circ\text{C}$ ):

PE DNA Thermal Cycler 480:

- 94°C for 1 min
- 25–30 cycles:
  - 94°C 30 sec
  - 68°C 4 min<sup>‡</sup>

PE GeneAmp Systems 2400/9600  
(or hot-lid thermal cycler):

- 94°C for 30 sec
- 25–30 cycles:
  - 94°C 5 sec
  - 68°C 4 min<sup>‡</sup>

**Notes on cycling:**

You may need to determine the optimal cycling parameters for your gene empirically.

<sup>‡</sup> The optimal extension time depends on the length of the fragment being amplified. We typically use 4 min for cDNA fragments of 2–5 kb. For 0.2–2-kb targets, we reduce the extension time to 2–3 min. For 5–10-kb targets, we increase the extension time up to 10 min.

5. When cycling is completed, analyze 5  $\mu\text{l}$  from each tube, along with appropriate DNA size markers, on a 1.2% agarose/EtBr gel.
6. [Optional] If the primary PCR reaction fails to give the distinct band(s) of interest or produces a smear, you may wish to perform a Southern blot using a cDNA probe or a nested primer as a probe.

Or, you may wish to perform a secondary, or “nested,” PCR reaction using the AP2 primer supplied with Marathon-Ready cDNA and a NGSP. (See the discussion on Primer Design in Section V of the User Manual.)

- a. Dilute 5  $\mu\text{l}$  of the primary PCR product into 245  $\mu\text{l}$  of Tricine-EDTA buffer.
- b. Repeat steps 1–5 above, using:
  - 5  $\mu\text{l}$  of the diluted primary PCR product in place of the Marathon-Ready cDNA.
  - 1  $\mu\text{l}$  of the AP2 primer and 1  $\mu\text{l}$  of your nested antisense GSP.
  - Fewer cycles (15–20 instead of 25–30).

## VII. Characterization of RACE Products

At this point, we recommend that you characterize your RACE fragments and confirm that you have amplified the desired product. This can be done by (1) comparison of RACE products obtained with GSPs and NGSPs; (2) Southern blotting; and (3) cloning and sequencing. Characterization of your RACE products at this point can prevent confusion and wasted effort in your subsequent experiments, even when both RACE reactions produce a single major product.

After RACE products have been characterized by partial or complete sequencing, you have two options for generating the full-length cDNA:

- 1) Generation of Full-Length cDNA by PCR (Section VIII).
- 2) Generation of Full-Length cDNA by Cloning (Section IX).

**Notice to Purchaser**

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes, nor is it intended for human use. Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Clontech Laboratories, Inc.

Suppression PCR is covered by U.S. Patent No. 5,565,340. Foreign patents pending.

A license under U.S. Patent Nos. 4,683,202, 4,683,195, and 4,965,188 and U.S. Patents Nos. 5,407,800, 5,322,770, and 5,310,652 or their foreign counterparts, owned by Roche Molecular Systems, Inc, and F. Hoffmann-La Roche Ltd, for use in research and development, has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction (PCR) and related processes described in said patents where the processes are covered by patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end-user in order to have a complete license to use this product in the PCR process where the process is covered by patents. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an Authorized Thermal Cycler. No right to perform or offer commercial services of any kind using PCR, where the process is covered by patents, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR Process where the process is covered by patents may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

Clontech, Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.  
Clontech is a Takara Bio Company. ©2005